

## SERVICE PACK CEREC SW 4.3.1

# Update Description

This service pack CEREC SW 4.3.1 updates the existing CEREC SW 4.3.0 installation. It supports both CEREC Omnicam and CEREC Bluecam and contains the following the following changes and corrections:

### Changes

- The \*.ssi data set file size has been reduced permitting easier importation into Sirona GALAXIS software.
- Omnicam models in 3D PDF format are now shown in color. If available, the SmileDesign face in 3D PDF can also be displayed.
- The base of the restoration model is now shown in a lighter shade to increase the contrast to the model.
- The color gradient can be regulated in Menu Analysis Tools and also deactivated.
- The color gradient for certain restoration materials has been activated. Formerly grey in color.
- Block sizes that are too small for the designed restoration will no longer be hidden in the step “Change Block Size” but now also indicated.
- Default has been changed for export of data. (YYYY-MM-DD\_indication dental notation.xyz; e.g. 2014-08-07\_I45\_C36\_V43\_C11\_C22\_C25\_I16\_C24\_FDI.rst)
- Sirona inCoris TZI C has been included as a new material.
- Ivoclar TelioCAD meso has been included as a new material for abutments
- Ivoclar Empress CAD & Empress CAD multi have been activated as a veneering structure for abutments.

### Corrections

- Errors occurring under certain circumstances when using Incisal Variation tools have been rectified.
- The problem with artifacts appearing when using the Shape tools has been corrected.
- The effect occurring under certain circumstances “Hide Neighbors“ where the neighbor was not available has been fixed.
- Minor bugs in the virtual articulator have been fixed and the calculation speed increased.
- An error in the Omnicam calibration process has been rectified.
- The problem that the Connect Button was not always available has been corrected.
- The preparation line for pontics can once again be edited through the transparent Biocopy Model.
- In Settling the icon images of the buttons Erase & Add show now the correct functionality.
- One can now sort in the Patient List.

We recommend all CEREC SW 4.3 users to install the service pack.