# Sidexis 4

- Installationsvoraussetzungen V4.3.1
- Installation Requirements V4.3.1
- Conditions d' installation V4.3.1
- Requisitos de instalación V4.3.1
- Requisiti di installazione V4.3.1

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# Dear customer,

We are delighted that you chose to place your trust in Dentsply Sirona's innovations and wish to congratulate you on the purchase of your new dental product.

Your product is marked with the adjacent symbol. This means that the instructions for use are made available in electronic form on our online portal.

You can easily download the Operating Instructions online from <http://www.dentsplysirona.com/ifu>. If you use a terminal device which functions with a QR app, you have the additional advantage of calling up the Operating Instructions anytime, anywhere. Simply scan the QR code and it will connect you directly to our online portal.

If you would prefer to have the instructions for use in paper format, please request them at <http://www.dentsplysirona.com/ifu>. Simply fill out the form online and we will be delighted to send you a printed copy, free of charge. The instructions for use will be received within the next 7 days. We would like to thank you in advance for your cooperation.

We hope that you thoroughly enjoy using your Dentsply Sirona products.

Your Dentsply Sirona Team

<http://www.dentsplysirona.com/ifu>

# Notes

This document covers the most important conditions that a practice must meet for the installation of Sidexis 4 Version 4.3.1.

## General requirements for the installation of Sidexis 4

The current patch version is required for each operating system.

Note: Please keep your operating system and other Microsoft components up-to-date with the latest security updates, even after installation.

You must be connected to the Internet during the installation. For the installation, a Windows user with administrator rights must be logged in.

The system requirements described below apply to Sidexis 4 including the operating system and the accompanying components (such as MS SQL Server Express). If other software packages are run in parallel on the same computers, their system requirements must be added.

A LAN bandwidth of 1 GB/s is required for smooth operation on the network. Please note that parallel-running applications and services can reduce the bandwidth available for Sidexis 4 and thus affect the speed of data processing in the network.

Sidexis 4 requires a MS SQL Server. Sidexis 4 is compatible with MS SQL Server versions 2014 and 2017. The MS SQL Server 2014 Express (SP3) is included and installed by the Sidexis 4 setup, provided a compatible version is not already installed.

## Notes on virtualization:

Sidexis 4 is released for certain Windows operating systems. It is assumed that the released operating system runs in a stable and high-performance system environment. Whether this system environment is based on a virtual machine or a physical computer is the responsibility of the respective IT manager of the customer.

Services can be easily virtualized. The virtualization of graphics plays the most important role in virtualizing applications. The virtualization must provide smooth graphics performance when the maximum number of stations are running.

It must be noted that virtualization results in higher resource consumption. Therefore, setting only the minimum system requirements does not lead to a satisfactory operating behavior of applications. For this reason, the resources made available to the individual virtual machines must be delimited; we advise against dynamic settings. Sidexis 4 checks the avail-

able memory, for example, so that transactions may be aborted before the operating system can "dynamically" expand the available memory.

The performance thresholds of the virtualization used must be known to the operator and must not be exceeded in real operation.

## Hardware requirements

To meet the different requirements depending on the size of the practice or clinic, the recommended hardware requirements for the server PC are differentiated according to the number of workstations. This allows a better working speed (performance) to be achieved depending on the size of the practice.

The hardware requirements for the workstations are independent of the size of the practice. For information on how to achieve better performance on workstations, see the "Tips for better performance" section.

## Server PC

	Minimum requirements	recommended for up to 20 work stations	recommended for up to 50 work stations	recommended for more than 50 work stations
<b>Operating system</b>	Windows 8.1 Pro (64 bit) <sup>1</sup> Windows 10 Pro (64 bit) <sup>1,2</sup> Windows 10 Pro (64 bit) <sup>1,2</sup> Windows Server 2012 R2 Windows Server 2016 Windows Server 2019	Windows Server 2012 R2 Windows Server 2016 Windows Server 2019	Windows Server 2016 Windows Server 2019	Windows Server 2016 Windows Server 2019
<b>MS SQL server</b>	MS SQL Server 2014 Express (included in the Sidexis 4 setup)	MS SQL Server 2014 Express (included in the Sidexis 4 setup) or higher version	MS SQL Server 2014 or 2017 (not included in the scope of delivery)	MS SQL Server 2014 or 2017 (not included in the scope of delivery)
<b>CPU</b>	2 GHz / 2-core processor	2.3 GHz / 4-core processor	2.3 GHz / 4-core processor	2.4 GHz / 12-core processor (e.g. XEON E5-2673v3)
<b>RAM</b>	4 GB	8 GB	16 GB	32 GB
<b>Hard disk</b>	> 1 TB	> 1 TB	> 1 TB free space	> 1 TB
	Sufficient disk space must be ensured even during operation.			
<b>Graphics card</b>	any	any	any	any
<b>Network</b>	1 Gbps Ethernet with DNS	1 Gbps Ethernet with DNS	1 Gbps Ethernet with DNS	5-10 Gbps Ethernet with DNS

- Using these systems for server operation is not recommended, as they are not designed by Microsoft for this purpose. Slower performance compared to a server operating system is to be expected.
- Tested: Version 1903; Build 18362.175.

# Workstation PC

	Minimum for 2D stations	Minimum for 3D stations
<b>Operating system<sup>1</sup></b>	Windows 8.1 Pro (64 bit) Windows 10 Pro (64 bit) <sup>2</sup> Windows 10 Enterprise (64 bit)	
<b>CPU</b>	2 GHz / 2-core processor	2.3 GHz / 4-core processor
<b>RAM</b>	4 GB	8 GB
<b>Hard disk</b>	5 GB	5 GB
<b>Graphics card</b>	512 MB, DirectX 9.0c compatible or Intel onboard graphics <sup>3</sup>	2 GB, DirectX 10 compatible or Intel onboard graphics <sup>3</sup>
<b>Network</b>	1 Gbps Ethernet with DNS	
<b>Screen resolution</b>	1280x1024 pixels	

1. For the 64-bit operating systems, installation under bootcamp is also approved.
2. Tested: Version 1903; Build 18362.175.
3. To ensure that the interaction with the 3D rendered volume is reliably "judger-free", graphics cards are recommended with at least the following Passmark GPU benchmark values:  
NVIDIA: Passmark > 1000  
AMD: Passmark > 1500  
Intel Onboard: Passmark > 540

The Passmark values can be found at [www.videocardbenchmark.net](http://www.videocardbenchmark.net), provided the name of the graphics card is known.

## Tips for better performance

- More working memory (RAM) improves the performance for most tasks. RAM should be increased especially if other applications are run at the same time as Sidexis 4.
- Use an SSD instead of a conventional hard disk
- A good graphics card with its own graphics memory is advantageous for working with 3D data. The Intel on-board graphics card does not have its own graphics memory and uses the RAM of the PC. A corresponding increase of the RAM is therefore recommended.

# Compatibility

## Sidexis 4 practice devices

All devices in the current imaging systems portfolio are supported

Supported	No longer supported
Xios, Xios Plus, Xios XG, Xios AE	Sidexis IO sensors
XIOS Scan	Orthophos DS
Orthophos XG	Pan Elite
Orthophos XG 3D	
Orthophos SL, Orthophos S, Orthophos E	
Axeos	
Galileos	
FaceScan	

- Please observe the individual product installation requirements as well as the document "Supplements for operation with Sidexis 4" (download via dealer area, REF 65 24 768)
- Sidexis 4 is only tested and released with the respective current PC software components. An update of all PC software components to the latest version is required (may require an update of the device software, and thus recalibration).

## Cameras

Integration is possible for the cameras listed.

- SiroCam 3, C, digital, AF, AF+, F
- SiroCam UAF Plus

## Plugins

Plugins need to be installed to operate the cameras

- Video module: Available in the Sidexis 4 installation menu under "Accessories"
- Newer versions may be available for download in the dealer area under the menu item "Product Info -> Treatment Centers"

The individual product installation requirements must be met.

## Integration of Dentsply Sirona treatment centers:

Dental Unit Plugins are available for Sidexis 4. When connected to Dentsply Sirona treatment centers, they enable

- the connection of the Intraoral camera
- the additional display of the ApexLocator on a monitor\*
- the automatic time synchronization\*
- the patient-specific chair positions\*
- the automatic head rest adjustment\*

The Dental Unit Plugins can be found in the Sidexis 4 installation menu under "Accessories".

Newer versions may be available for download in the dealer area under the menu item

"Product Info -> Treatment Centers"

\*Not available for all Dentsply Sirona treatment centers

## Other devices

- VistaScan image plate scanners can be used via a plugin from Dürr Dental (Dürr Bridge).
- Third-party intraoral cameras can be connected to Sidexis 4 (version 4.1.3 or higher) via WDM (Windows Driver Model). To do this, install the video module from the Sidexis 4 setup.
- As of CEREC Version 4.4, the CEREC database can be connected to Sidexis 4 via a new interface.
- The Sidexis SyncServices of REINKE SOFTWARE ENGINEERING (RSWE) and Sidexis 4 are completely compatible as of version 4.1.3. This requires a current version of Sidexis SyncServices (available from RSWE).

## Sidexis extensions that are not yet supported

The following Sidexis XG extensions are not yet supported by Sidexis 4 version 4.3.1:

- TSE version

Should the customer continue to need these functions, then Sidexis XG is the more suitable product for these customers.

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Sirona Dental Systems GmbH