

CEREC SW 4.6.1.152739 & CEREC CAM SW 4.6.1

MATERIAL PACK DESCRIPTION

The free of charge material pack requires a pre-installed version CEREC SW 4.6.1.152739 for CEREC SW update and a pre-installed version CEREC CAM SW 4.6.1152708 for CEREC Cam SW update.

You need to install the free of charge update CEREC SW 4.6.1.152739 first, if your system currently runs with SW version CEREC SW 4.6.0 or a previous 4.6.1 SW version (CEREC SW 4.6.1.141333 und CEREC SW 4.6.1.148853). You find the file download for CEREC SW 4.6.1.152739 and the update description on my.cerec.com.

You need to install the free of charge update CEREC CAM SW 4.6.1. first, if your system currently runs with SW version CEREC CAM SW 4.6.0. You find the file download for CEREC CAM SW 4.6.1. on my.cerec.com.

News and changes with CEREC SW 4.6.1/CEREC CAM 4.6.1 with material pack:

- Integration of new materials:

Manufacturer	Material
GC	LiSi Block, Grinding & Sintering, Block 14, not for USA
Ivoclar	IPS e.max ZirCAD MT Multi, Grinding & Milling & Sintering, Block C17, Block B45, not worldwide available
Kuraray Noritake Dental	Avencia Block, worldwide available now; KATANA Zirconia Block, Block STML 14 Z L, not for USA, not for JPN
VITA	ENAMIC multiColor, Block EMC-16
3M	Chairside Zirconia, Block Crown, Block Bridge, Sintering; not worldwide available

- Optimized sintering process for KND KATANA Zirconia blocks.

Please observe regional restrictions regarding TiBase and materials.

Refer to document CEREC-4_6- TiBase-Material-Compatibility_EN.pdf at my.cerec.com.

Installation & License

- Automatic update:

A pop-up message informing about the new update shows up. Start the update process directly by pressing "Install" or start the process later.

When you like to start the update process later you have to go to the Windows taskbar. Open the „AutoUpdate@ Sirona Connect Center“ by clicking on the corresponding icon.



The CEREC SW 4.6.1 Material Pack / CEREC CAM SW 4.6.1 Material Pack is listed in the update list of the „AutoUpdate @ Sirona Connect Center“.

Start the download by clicking the button. If the network is disconnected the download can be continued at the same position by clicking the button again. The second button starts the installation setup.

- You have to unpack the update files first, if you download them from my.cerec.com. Afterwards you can start the unpacked .exe file in the main directory.
- No additional license needed to update from a pre-installed CEREC SW 4.6.x to CEREC SW 4.6.1.152739 with material pack.
- No additional license needed to update from a pre-installed CEREC CAM SW 4.6.x to CEREC CAM SW 4.6.1. with material pack.
- The SW version number 4.6.1 is kept. The successful installation of the material pack can be identified by looking at the included materials or by checking the block database version and process data base version in the extended about dialogue.



CEREC 4.6.1:

Block DB: 461.1.437 → 461.1.478; Process DB: 461.1.38 → 461.1.40

CEREC CAM 4.6.1:

Block DB: 461.64.444 → 461.64.477; Process DB: 461.64.38 → 461.64.40