

# Reservoir Waterline Maintenance Steps

## Cavitron Select® SPS™ Unit with Cavitron Select® Reservoir

Follow these steps to maintain the waterline of a Cavitron Select® SPS™ with Cavitron Select® Reservoir System.

1



After completing the “Shut-Down Procedures (end of day)” referenced in your Cavitron Select Instructions for Use, begin by preparing a 1:10 Sodium Hypochlorite Solution by mixing 1 part 5.25% Sodium Hypochlorite with 10 parts tap water.

Remove reservoir cap and pour the 1:10 Sodium Hypochlorite Solution into the Cavitron Select® SPS™ Reservoir Bottle. **Do not overfill with solution.**

2



Securely replace reservoir cap on bottle.

3



Turn the System ON.

Adjust the lavage/water control at the base of the handpiece to the maximum water flow.


4



Begin to flush the 1:10 Sodium Hypochlorite Solution through the system for at least 30 seconds.



5



Turn the System OFF. The system should be left undisturbed, allowing the prepared solution to sit in the waterline for a **minimum of 10 minutes but no more than 30 minutes.**

It is recommended that a sign be placed on the scaling system: “Do Not Use System is Being Treated”





After the minimum **10-minute wait time**, unscrew reservoir cap and discard disinfectant solution.



Fill Cavitron Select® SPS™ Reservoir Bottle with tap or distilled water. **Do not overfill with liquid.**



Turn the System ON.

Flush clean water through the system for **2 minutes** or until all Sodium Hypochlorite odor disappears.



Turn the system OFF.

Detach Cavitron Select® SPS™ scaler from Select Reservoir System.



Detach Cavitron® Select Reservoir Bottle from Pump.



Autoclave Cavitron Select SPS™ Reservoir Bottle per Infection Control Procedure booklet supplied with your system.

Ensure that all the Sodium Hypochlorite solution has been completely flushed from the system prior to use on a patient.

**Your system is now ready for use.** Follow the CDC's guidelines for the frequency of testing and perform testing to your Cavitron Select® SPS™ Scaling System.