

Achieve your best smile
with the SureSmile® Aligner
system, now including
SureSmile® VPro™

1

Consult with your doctor
and begin a SureSmile®
Clear Aligner programme

2

Use your SureSmile VPro
device for 5 minutes per day

3

Smile big and smile often

Download the SureSmile®
VPro™ App

Pair the SureSmile® VPro™ app with your
SureSmile® VPro™ device in order to:

- Track your daily SureSmile® VPro™ usage
- Receive reminders to use device daily
and change aligners as prescribed by
your doctor
- Monitor SureSmile® VPro™ battery life
- Share results and progress photos



Ask your doctor about
SureSmile® Clear Aligners
and SureSmile® VPro™



Learn more at [SureSmile.com](https://www.suresmile.com)

© 2023 Dentsply Sirona. All rights reserved.

Intended for use by the orthodontic patient during treatment with aligners.
Rx only. See Instructions for Use for full prescribing information.



Dentsply Sirona, PO Box 2024, SE-128 30 SKARPNÄCK
www.dentsplysirona.com/suresmile
All Rights Reserved.



The little device
with a big impact



SureSmile® VPro™
High-Frequency Vibration Device



Help meet your treatment goals

SureSmile® VPro™ delivers gentle vibrations that help aligners fit properly.*

When used with your aligners for 5 minutes a day, your SureSmile VPro Device:

- May reduce discomfort^{1,2}
- May reduce office visits¹
- May speed up treatment time by up to 50%^{1,3}

Intended for use by the orthodontic patient during treatment with aligners. Rx only. See Instructions for Use for full Prescribing Information.

*Over not using the device.

References: 1. Alansari S, Atique MI, Gomez JP, et al. The effects of brief daily vibration on clear aligner orthodontic treatment, *J World Fed Ortho*. 2018. <https://doi.org/10.1016/j.ejwf.2018.10.002>. 2. Shipley T. Effects of high frequency acceleration device on aligner treatment-A Pilot Study. *Dent J*. 2018;6(3):32. 3. Alikhani M, Alansari S, Hamidaddin MA, Sangsuwon C, Alyami B, et al. Vibration paradox in orthodontics: anabolic and catabolic effects. *PLoS ONE*. 2018;13(5):e01965.



Results worth smiling about

CASE EXAMPLES



Estimated treatment time: 19 months*

Actual treatment time: 9 months



Estimated treatment time: 19 months*

Actual treatment time: 6 months



Estimated treatment time: 12 to 14 months*

Actual treatment time: 5 months

Patient results may vary. Case examples courtesy of Dr David Boschken, Dr Bill Dischinger, and Dr David Boschken, respectively. Doctors are located in the US and may have been compensated for their participation.

*Estimated treatment times provided by treating clinician and calculated with 2-week aligner intervals. Cases may have been treated with an earlier VPro model featuring a variety of clear aligner systems. Case examples are from real patients. Individual results may vary and are not intended to represent or guarantee that anyone will achieve the same or similar results. See your treating dental professional for advice unique to you and your oral health.