

SureSmile® VPro™ Quick Start Guide



Why SureSmile VPro?

SureSmile VPro delivers gentle vibrations that help aligners fit properly.[†]

When used with your aligners for 5 minutes a day, your SureSmile VPro device:

- **May reduce discomfort**^{1,2}
- **May reduce office visits**¹
- **May speed up treatment time by up to 50%**^{1,3}

How to Start?



Rinse the mouthpiece under water prior to first use
Make sure the mouthpiece is detached from vibration device before rinsing



Charge your device fully before your first session



Begin 5-minute session and wear aligners during use

Intended for use by the orthodontic patient during treatment with aligners to facilitate minor anterior tooth movement. Rx only. See Instructions for Use for full prescribing information.

[†]Over not using the device.

References: 1. Alansari S, Atique MI, Gomez JP, et al. The effects of brief daily vibration on clear aligner orthodontic treatment, *J World Fed Ortho*. 2018. <https://doi.org/10.1016/j.ejwof.2018.10.002>. 2. Shipley T. Effects of high frequency acceleration device on aligner treatment-A Pilot Study. *Dent J*. 2018;6(3):32. 3. Alikhani M, Alansari S, Hamidaddin MA, Sangsuwon C, Alyami B, et al. Vibration paradox in orthodontics: anabolic and catabolic effects. *PLoS ONE*. 2018;13(5):e01965.

SureSmile® VPro™ Quick Start Guide

Light Notifications

Blue (3 blinks): Session starts

Blue (slow blinking): Charging in progress

Green (1 blink): App has connected

Green (3 blinks): Session ends

Green (slow blinking): Charging complete

Green (blinks and vibrates): Device locate

Yellow (2 blinks): Successful app pairing

Red (slow blinking): Low battery

Red, Green, Blue (1 blink each): Device reset

Frequently Asked Questions

■ How much pressure should I apply?

Bite down gently on the mouthpiece during use. There is no need to bite down harder than as needed to support the device.

■ How do I clean my mouthpiece?

Remove the mouthpiece from the vibration device and rinse the mouthpiece under water, brushing with a toothbrush as necessary. Keep the vibration device dry.

■ What if my SureSmile VPro won't pair or sync with the app?

First, make sure that you have your Bluetooth® connection turned on. Next, make sure the device is charged and in close proximity to your mobile device. Close and reopen the app. If you are still experiencing difficulty, press and hold the button on the device for 10 seconds to reset the device.

Download the SureSmile® VPro™ App

Pair the SureSmile VPro app with your SureSmile VPro device in order to:

- Track your daily SureSmile VPro usage
- Receive reminders to use device daily and change aligners as prescribed by your doctor
- Monitor your SureSmile VPro battery life
- Share results and progress photos



Learn more at [SureSmile.com](https://www.suresmile.com)

©2023 Dentsply Sirona. All Rights Reserved. 202301. APAC/CAN/EU/INDIA/US