

# inPost compatibility overview for scanning on multi-unit abutments from Medentika

Implant manufacturer	Implant system	Platform	Multi-Unit Abutments	Medentika: MedentiBASE (Gingival height in mm)					SIRONA	
			Series	0,5	1,5	2,5	3,5	4,5	inPost	REF.
<b>BioHorizons</b>										
	Tapered plus, Tapered internal, Tapered internal tissue level, Internal dental implant, Single stage dental implants	<span style="color: yellow;">■</span> 3.5 <span style="color: green;">■</span> 4.5	R-Series	R4000 R4010	R4100 R4110	R4200 R4210	R4300 R4310	R4400 R4410	inPost for MedentiBASE Standard inPost for MedentiBASE Standard	6552298 6552298
<b>Biomet 3i</b>										
	NanoTite Parallel Walled, OSSEOTITE Parallel Walled, OSSEOTITE XP, Full OSSEOTITE Parallel Walled, Full OSSEOTITE XP, NanoTite Tapered, OSSEOTITE Tapered, Full OSSEOTITE Tapered	<span style="color: purple;">■</span> 3.4 <span style="color: blue;">■</span> 4.1	I-Series			I4000 I4010	I4100 I4110	I4200 I4210	inPost for MedentiBASE Standard inPost for MedentiBASE Standard	6552298 6552298
	NanoTite Certain® PREVAIL, OSSEOTITE Certain PREVAIL, NanoTite Tapered Certain PREVAIL, NanoTite Parallel Walled Certain, OSSEOTITE Parallel Walled Certain, Full OSSEOTITE Parallel Walled Certain, OSSEOTITE XP Certain, Full OSSEOTITE XP Certain, NanoTite Tapered Certain, OSSEOTITE Tapered Certain, Full OSSEOTITE Tapered Certain	<span style="color: purple;">■</span> 3.4 <span style="color: blue;">■</span> 4.1	H-Series	H4000 H4010	H4100 H4110	H4200 H4210	H4300 H4310	H4400 H4410	inPost for MedentiBASE Standard inPost for MedentiBASE Standard	6552298 6552298
<b>Bredent</b>										
	Sky	<span style="color: yellow;">■</span> narrow 3.5	B-Series	B4000	B4100	B4200	B4300	B4400	inPost for MedentiBASE Standard	6552298
<b>Camlog</b>										
	Camlog	<span style="color: grey;">■</span> 3.3 <span style="color: yellow;">■</span> 3.8 <span style="color: red;">■</span> 4.3	C-Series	C4000 C4010 C4020	C4100 C4110 C4120	C4200 C4210 C4220	C4300 C4310 C4320	C4400 C4410 C4420	inPost for MedentiBASE Standard inPost for MedentiBASE Standard inPost for MedentiBASE Standard	6552298 6552298 6552298
	Conelog	<span style="color: grey;">■</span> 3.3 <span style="color: yellow;">■</span> 3.8	D-Series	D4000 D4010	D4100 D4110	D4200 D4210	D4300 D4310	D4400 D4410	inPost for MedentiBASE Standard inPost for MedentiBASE Standard	6552298 6552298
<b>Dentsply Implants (AstraTech, Friadent)</b>										
Astra Tech	OsseoSpeed TX	<span style="color: teal;">■</span> 3.5/4.0 <span style="color: purple;">■</span> 4.5/5.0	S-Series	S4000 S4020	S4100 S4120	S4200 S4220	S4300 S4320	S4400 S4420	inPost for MedentiBASE Standard inPost for MedentiBASE Standard	6552298 6552298
Friadent	Frialit / Xive	<span style="color: grey;">■</span> 3.4 <span style="color: yellow;">■</span> 3.8 <span style="color: blue;">■</span> 4.5	T-Series	T4000 T4005 T4010	T4100 T4105 T4110	T4200 T4205 T4210	T4300 T4305 T4310	T4400 T4405 T4410	inPost for MedentiBASE Standard inPost for MedentiBASE Standard inPost for MedentiBASE Standard	6552298 6552298 6552298
	Ankylos	3.5-7.0	Y-Series	Y4000	Y4100	Y4200	Y4300	Y4400	inPost for MedentiBASE Standard	6552298
<b>MIS</b>										
	Seven (internal hex.)	<span style="color: purple;">■</span> Standard <span style="color: green;">■</span> Wide	R-Series	R4000 R4010	R4100 R4110	R4200 R4210	R4300 R4310	R4400 R4410	inPost for MedentiBASE Standard inPost for MedentiBASE Standard	6552298 6552298
<b>Nobel Biocare</b>										
	<b>NobelReplace®</b> (Straight, Tapered, Platform Switching), <b>Replace Select</b> (Straight, Tapered, TC), <b>NobelSpeedy Replace</b>	<span style="color: red;">■</span> NP <span style="color: yellow;">■</span> RP	E-Series	E4000 E4010	E4100 E4110	E4200 E4210	E4300 E4310	E4400 E4410	inPost for MedentiBASE Standard inPost for MedentiBASE Standard	6552298 6552298
	Nobel Active, NobelReplace Conical Connection	<span style="color: red;">■</span> NP <span style="color: yellow;">■</span> RP	F-Series	F4000 F4010	F4100 F4110	F4200 F4210	F4300 F4310	F4400 F4410	inPost for MedentiBASE Standard inPost for MedentiBASE Standard	6552298 6552298
	Brånemark® Mk III, Mk IV NobelSpeedy Groovy	NP RP	K-Series			K4000 k4010	K4100 K4110	K4200 K4210	inPost for MedentiBASE Standard inPost for MedentiBASE Standard	6552298 6552298
<b>Straumann</b>										
	Standard (Roxolid), Standard Plus, Tapered Effect	NNC (3,5 mm) RN (4,8 mm)	N-Series	N4030 N4010	N4130	N4230	N4330	N4430	inPost for MedentiBASE Standard inPost for MedentiBASE N4010	6552298 6561091
	Bone Level	NC (3.3 mm) RC (4,1/4,8 mm)	L-Series	L4000 L4010	L4100 L4110	L4200 L4210	L4300 L4310	L4400 L4430	inPost for MedentiBASE Standard inPost for MedentiBASE Standard	6552298 6552298
<b>Zimmer</b>										
	Tapered Screw-Vent	<span style="color: green;">■</span> 3.5 <span style="color: purple;">■</span> 4.5	R-Series	R4000 R4010	R4100 R4110	R4200 R4210	R4300 R4310	R4400 R4410	inPost for MedentiBASE Standard inPost for MedentiBASE Standard	6552298 6552298

Technische Änderungen und Irrtümer vorbehalten, Best.-Nr. A91100-B944-01-7600, Printed in Germany, Dispo-Nr. 04605\_1015.V0