


CEREC SW 5.2.4

CEREC MC XL Premium Package, CEREC MC XL, CEREC MC X

| Motor set | Left | Right | Combination |
|------------------|--------------------|----------------------------|------------------------------|
| Set1 (or Set2**) | Step Bur 12S | Cylinder Pointed Bur 12S | 1 |
| Set1 (or Set2**) | Step Bur 12 | Cylinder Pointed Bur 12S | 2 |
| Set1 (or Set2**) | Step Bur 20 | Cylinder Pointed Bur 20 | 3 |
| Set1 | Shaper 25 RZ | Finisher 10 | 4 ⁽¹⁾ Wet Milling |
| Set1 | Cylinder Bur 12 EF | Cylinder Pointed Bur 12 EF | 5 ⁽²⁾ Extra Fine |
| Set1 | Shaper 25 | Finisher 10 | 6 ⁽³⁾ Dry Milling |

⁽¹⁾ Only with grinding/milling units from the following serial numbers: inLab MC XL 129001, CEREC MC XL 129001, CEREC MC XL Premium Package 302001, CEREC MC X 231001 or replacement of the left motor on 2-motor machines or of the left motor in the motor set 1 on 4-motor machines.
⁽²⁾ Only for 4-motor grinding/milling units
⁽³⁾ Requires CEREC grinding/milling unit with dry milling configuration and CEREC suction unit.
⁽⁴⁾ USA/Japan
⁽⁵⁾ Worldwide except Japan
⁽⁶⁾ Japan
⁽⁷⁾ Worldwide except USA
⁽⁸⁾ Not available in all countries

| Description | Figure | Ref. |
|----------------------------|--------|-----------|
| Step Bur 12S | 1 | 62 40 167 |
| Step Bur 12 | 2 | 62 60 025 |
| Cylinder Pointed Bur 12S | 3 | 62 40 159 |
| Cylinder Bur 12 EF | 4 | 65 35 186 |
| Cylinder Pointed Bur 12 EF | 5 | 65 35 178 |
| Step Bur 20 | 6 | 62 59 597 |
| Cylinder Pointed Bur 20 | 7 | 62 59 589 |
| Shaper 25 RZ | 8 | 64 33 440 |
| Finisher 10 | 9 | 62 99 387 |
| Shaper 25 | 10 | 62 99 395 |



| Manufacturer | Material | Combination |
|------------------|--|-------------|
| Dentsply Sirona | CEREC Tessera | 1 2 5 |
| | CEREC Tessera Abutment | 1 2 5 |
| | CELTRA® Duo | 1 2 5 |
| | CEREC Blocs C / CEREC Blocs | 1 2 5 |
| | CEREC Blocs C PC / CEREC Blocs PC | 1 2 5 |
| | CEREC Blocs C In | 1 2 5 |
| | CEREC Guide Blocs ⁽⁷⁾ (for CEREC Guide 2 und 3) | 4 |
| | CEREC Zirconia | 3 4 6 |
| | CEREC Zirconia meso | 3 4 6 |
| | CEREC Zirconia+ | 4 6 |
| | CEREC MTL Zirconia | 4 6 |
| | inCoris TZI C | 3 4 6 |
| | inCoris TZI | 3 4 6 |
| | inCoris ZI meso | 3 4 6 |
| 3M | Lava™ Ultimate | 1 2 5 |
| | Paradigm™ MZ 100 ⁽⁴⁾ | 1 2 5 |
| | Chairside Zirconia | 6 |
| COLTENE | BRILLIANT Crios | 1 2 5 |
| GC | CERASMART™ 270 | 1 2 5 |
| | CERASMART™ 300 ⁽⁶⁾ | 1 2 5 |
| | CERASMART™ LAYER ⁽⁶⁾ | 1 2 5 |
| | CERASMART™ PRIME ⁽⁶⁾ | 1 2 5 |
| | Initial™ LRF | 1 5 |
| | Initial™ LiSi Block ⁽⁷⁾ | 1 2 |
| MERZ DENTAL | artBloc® Temp | 1 4 |
| KURARAY NORITAKE | KATANA™ Zr STML | 3 4 6 |
| | KATANA™ Avencia Block | 1 2 5 |
| | KATANA™ Avencia N ⁽⁶⁾ | 1 2 5 |

| Manufacturer | Material | Combination |
|----------------------------------|------------------------------------|-------------|
| SHOFU | BLOCK HC | 1 2 5 |
| | BLOCK HC (2L) | 1 2 5 |
| | BLOCK HC HARD ⁽⁸⁾ | 1 2 5 |
| | BLOCK HC HARD AN ⁽⁶⁾ | 1 2 5 |
| | BLOCK HC SUPERHARD ⁽⁶⁾ | 1 2 5 |
| | BLOCK ZR LUCENT ⁽⁶⁾ | 6 |
| VITA | VITABLOCS® MarkII | 1 2 5 |
| | VITABLOCS® TriLuxe | 1 2 5 |
| | VITABLOCS® TriLuxe forte | 1 2 5 |
| | VITABLOCS RealLife® | 1 2 5 |
| | VITA ENAMIC® | 1 2 5 |
| | VITA ENAMIC® IS | 1 5 |
| | VITA ENAMIC® multiColor | 1 2 5 |
| | VITA SUPRINITY® PC | 1 2 5 |
| | VITA CAD-Temp® monoColor | 1 4 |
| | VITA CAD-Temp® multiColor | 3 4 |
| VITA CAD-Temp® IS ⁽⁷⁾ | 3 4 | |
| VITA YZ HT | 3 4 6 | |
| IVOCLAR | IPS Empress® CAD | 1 2 5 |
| | IPS Empress® CAD Multi | 1 2 5 |
| | IPS e.max® CAD | 1 2 5 |
| | IPS e.max® ZirCAD LT | 3 4 6 |
| | IPS e.max® ZirCAD MT Multi | 3 4 6 |
| | IPS e.max® CAD Abutment | 1 5 |
| | Telio® CAD | 1 4 |
| | Telio® CAD Abutment ⁽⁷⁾ | 3 4 |
| | Tetric® CAD | 1 2 5 |

Order no. 62 38 369, Change no. 133197