



















Compatibility overview inPost for scanning on Multi Unit Abutments from nt-trading

			Multi-Unit Abutments	nt-trading: 2-CONnect	Dentsply Sirona	
Manufacturer	Implant system	Platform	Series		inPost	REF
BioHorizons	Tapered plus, Tapered internal, Tapered internal tissue level, Internal dental implant, Single stage dental implants	 3.5	R series	R 800 S R 810 S	inPost for 2-CONnect KS60	65 59 483
		 4.5				
Biomet 3i	NanoTite Certain® PREVAIL, OSSEOTITE Certain PREVAIL, NanoTite Tapered Certain PREVAIL, NanoTite Parallel Walled Certain, OSSEOTITE Parallel Walled Certain, Full OSSEOTITE Parallel Walled Certain, OSSEOTITE XP Certain, Full OSSEOTITE XP Certain, NanoTite Tapered Certain, OSSEOTITE Tapered Certain, Full OSSEOTITE Tapered Certain	 3.4	H series	H 800 S H 810 S	inPost for 2-CONnect KS60	65 59 483
		 4.1				
Dentsply Implants	Astra Tech: OsseoSpeed TX	 3.5/4.0  4.5/5.0	S series	S 800 S S 820 S	inPost for 2-CONnect KS61 inPost for 2-CONnect KS61	65 51 639 65 51 639
	Friadent: Frialit/Xive	 3.4  3.8  4.5	T- series	T 800 S T 805 S T 810 S	inPost for 2-CONnect KS60 inPost for 2-CONnect KS61 inPost for 2-CONnect KS61	65 59 483 65 51 639 65 51 639
MIS	Seven (internal hex.)	 Standard  Wide	R series	R 800 S R 810 S	inPost for 2-CONnect KS60 inPost for 2-CONnect KS60	65 59 483 65 59 483
Nobel Biocare	NobelReplace® (Straight, Tapered, Platform Switching), Replace Select (Straight, Tapered, TC), NobelSpeedy Replace	 RP  WP	E series	E 810 S E 820 S	inPost for 2-CONnect KS61 inPost for 2-CONnect KS61	65 51 639 65 51 639
	Nobel Active NobelReplace Conical Connection	 NP  RP	F series	F 800 S F 810 S	inPost for 2-CONnect KS61 inPost for 2-CONnect KS61	65 51 639 65 51 639
Osstem/Hiossen	Osstem TS (USA: Hiossen ET)	 Standard USA:Regular	S series	S 820 S	inPost for 2-CONnect KS61	65 51 639
Straumann	Standard (Roxolid), Standard Plus, Tapered Effect	RN (4,8 mm) WN (6,5 mm)	N series	N 810 S N 820 S	inPost for 2-CONnect KS61 inPost for 2-CONnect KS61	65 51 639 65 51 639
	Bone Level	NC (3,3 mm) RC (4,1/4,8 mm)	L series	L 800 S L 810 S	inPost for 2-CONnect KS61 inPost for 2-CONnect KS61	65 51 639 65 51 639
Zimmer	Tapered Screw-Vent	 3.5	R series	R 800 S R 810 S	inPost for 2-CONnect KS60	65 59 483
		 4.5				