


CEREC SW 5.2.3

CEREC MC XL Premium Package, CEREC MC XL, CEREC MC X

Motor set	Left	Right	Combination
Set1 (or Set2**)	Step Bur 12S	Cylinder Pointed Bur 12S	1
Set1 (or Set2**)	Step Bur 12	Cylinder Pointed Bur 12S	2
Set1 (or Set2**)	Step Bur 20	Cylinder Pointed Bur 20	3
Set1	Shaper 25 RZ	Finisher 10	4 ⁽¹⁾ Wet Milling
Set1	Cylinder Bur 12 EF	Cylinder Pointed Bur 12 EF	5 ⁽²⁾ Extra Fine
Set1	Shaper 25	Finisher 10	6 ⁽³⁾ Dry Milling

⁽¹⁾ Only with grinding/milling units from the following serial numbers: inLab MC XL 129001, CEREC MC XL 129001, CEREC MC XL Premium Package 302001, CEREC MC X 231001 or replacement of the left motor on 2-motor machines or of the left motor in the motor set 1 on 4-motor machines.
⁽²⁾ Only for 4-motor grinding/milling units
⁽³⁾ Requires CEREC grinding/milling unit with dry milling configuration and CEREC suction unit.
⁽⁴⁾ USA/Japan
⁽⁵⁾ Worldwide except Japan
⁽⁶⁾ Japan
⁽⁷⁾ Worldwide except USA
⁽⁸⁾ Not available in all countries

Description	Figure	Ref.
Step Bur 12S	1	62 40 167
Step Bur 12	2	62 60 025
Cylinder Pointed Bur 12S	3	62 40 159
Cylinder Bur 12 EF	4	65 35 186
Cylinder Pointed Bur 12 EF	5	65 35 178
Step Bur 20	6	62 59 597
Cylinder Pointed Bur 20	7	62 59 589
Shaper 25 RZ	8	64 33 440
Finisher 10	9	62 99 387
Shaper 25	10	62 99 395



Manufacturer	Material	Combination
Dentsply Sirona	CEREC Tessera	1 2 5
	CEREC Tessera Abutment	1 2 5
	CELTRA® Duo	1 2 5
	CEREC Blocs C / CEREC Blocs	1 2 5
	CEREC Blocs C PC / CEREC Blocs PC	1 2 5
	CEREC Blocs C In	1 2 5
	CEREC Guide Blocs ⁽⁷⁾ (for CEREC Guide 2 und 3)	4
	CEREC Zirconia	3 4 6
	CEREC Zirconia meso	3 4 6
	CEREC Zirconia+	4 6
	CEREC MTL Zirconia	4 6
	inCoris TZI C	3 4 6
	inCoris TZI	3 4 6
	inCoris ZI meso	3 4 6
3M	Lava™ Ultimate	1 2 5
	Paradigm™ MZ 100 ⁽⁴⁾	1 2 5
	Chairside Zirconia	6
COLTENE	BRILLIANT Crios	1 2 5
GC	CERASMART™ 270	1 2 5
	CERASMART™ 300 ⁽⁶⁾	1 2 5
	CERASMART™ LAYER ⁽⁶⁾	1 2 5
	CERASMART™ PRIME ⁽⁶⁾	1 2 5
	Initial™ LRF	1 5
	Initial™ LiSi Block ⁽⁷⁾	1 2
MERZ DENTAL	artBloc® Temp	1 4
KURARAY NORITAKE	KATANATM Zr STML	3 4 6
	KATANATM Avencia Block	1 2 5
	KATANATM Avencia N ⁽⁶⁾	1 2 5

Manufacturer	Material	Combination
SHOFU	BLOCK HC	1 2 5
	BLOCK HC (2L)	1 2 5
	BLOCK HC HARD ⁽⁸⁾	1 2 5
	BLOCK HC HARD AN ⁽⁶⁾	1 2 5
	BLOCK HC SUPERHARD ⁽⁶⁾	1 2 5
	BLOCK ZR LUCENT ⁽⁶⁾	6
VITA	VITABLOCS® MarkII	1 2 5
	VITABLOCS® TriLuxe	1 2 5
	VITABLOCS® TriLuxe forte	1 2 5
	VITABLOCS RealLife®	1 2 5
	VITA ENAMIC®	1 2 5
	VITA ENAMIC® IS	1 5
	VITA ENAMIC® multiColor	1 2 5
	VITA SUPRINITY® PC	1 2 5
	VITA CAD-Temp® monoColor	1 4
	VITA CAD-Temp® multiColor	3 4
VITA CAD-Temp® IS ⁽⁷⁾	3 4	
VITA YZ HT	3 4 6	
IVOCLAR VIVADENT	IPS Empress® CAD	1 2 5
	IPS Empress® CAD Multi	1 2 5
	IPS e.max® CAD	1 2 5
	IPS e.max® ZirCAD LT	3 4 6
	IPS e.max® ZirCAD MT Multi	3 4 6
	IPS e.max® CAD Abutment	1 5
	Telio® CAD	1 4
	Telio® CAD Abutment ⁽⁷⁾	3 4
	Tetric® CAD	1 2 5

Order no. 62 38 369; Change no. 132159