



# Protect your patients while polishing.

Are you complying with  
CDC Infection Prevention Guidelines  
when you polish?



The CDC Sterilization Guidelines state:

Dental handpieces and associated attachments, including low-speed motors and reusable prophylaxis angles, should always be heat sterilized between patients and not high-level or surface disinfected. <sup>1</sup>

The Nupro Freedom™ Cordless Prophy System with SmartMode™ Technology provides a cost-effective solution to comply with the CDC Infection Prevention Guidelines. Reduce your patients' risk of infection—and save your practice from the high costs of purchasing multiple handpieces.



Nupro Freedom® Disposable Prophy Angles are designed as a single-use disposable device to reduce the risk of contamination of patient material.

Disposa-shield® Disposable Polyethylene Barriers facilitate cleaning and protect the inner module from debris.

The CDC states that disposable prophy angles should be disposed of properly after use on one patient.

The FDA states to use a new polyethylene barrier for each patient and wipe down the Inner Module with a tuberculocidal disinfectant.

**SEALED COUPLING SHAFT:**  
The outer sheath is designed with a sealed coupling shaft which connects to the prophy angle on the outer sheath and a separate second coupling mechanism inside the sheath that connects to the inner module. This is uniquely different to the slip fit shaft design on other air driven low speed handpieces.



The Autoclavable Outer Sheaths have a fluid seal — isolating the inner module from patient contact and providing a secure infection control barrier.

The CDC states that dental handpieces should always be heat sterilized between patients.



Nupro® Prophy Paste comes in a disposable, unit dose cup to reduce the risk of patient cross-contamination. The Nupro family of prophy paste delivers a variety of options for patient preferences and needs.

Protect yourself and your patients with our Infection Prevention solutions.



Com-Fit Plus™ Natural Fit Masks are hypoallergenic face masks that offer comfort and breathability without sacrificing any protection capabilities. Our innovative elastic chin strip provides a better fit and helps improve protection.

We recommend using ASTM-2 for procedures with low to moderate fluid, aerosols and/or sprays, such as during the polishing procedure.



REVOLUTIONARY CHIN STRIP DESIGN



FlashTips® Disposable Air/Water Syringe Tips offer quick and easy changeover after each patient. With no converters required (on most common syringes), the disposable tips help reduce risk of cross-contamination.



CLOSE-UP OF INNER AND OUTER METAL TUBES



Signs of corrosion on inner tube

Your metal tip's tiny chamber could be home to nasty buildup from exposure to blood, saliva and dental materials. Over time, that buildup gets harder and harder to properly clean. Don't take the risk!

The CDC strongly recommends single-use disposable instruments whenever possible!



VoloWipes® Disinfecting Wipes have an alcohol and quaternary ammonium formula, which is an effective disinfectant, cleaner and deodorizer all in one.

Approved for cleaning and disinfecting Nupro Freedom's inner module, charging base, optional foot pedal, and power supply.



The CDC recommends wipe disinfection using a tuberculocidal surface disinfectant on the Nupro Freedom inner module. VoloWipes have been found effective in achieving this desired level of disinfection.

# Resources and support to help you deliver better, safer, faster infection prevention while polishing.

Instructions for Steam Sterilizing the Nupro Freedom Outer Sheath



Download here:  
[nuprofreedom.com](http://nuprofreedom.com)

Earn free CE credits.  
Pathways to Infection Prevention



Register here:  
[sultanacademy.com](http://sultanacademy.com)

Our online learning portal provides the knowledge, skills, and certification dental professionals need to stay up-to-date.



Find CE Courses, Technical Courses, and Practice Management Resources at:  
[dentsplysirona.com/academy](http://dentsplysirona.com/academy)

## Save money while complying with CDC Infection Prevention Guidelines.







### Nupro Freedom™ Cordless Prophy Package

- 1 Inner Module
- 3 Metal Outer Sheath
- 1 Charging Base
- 1 Handpiece Tray Cradle
- 1 Power Supply
- 1 Operation Manual
- 1 Box of Disposa-shield® Disposable Polyethylene Barriers (500 ct)
- 1 Box Nupro Freedom™ Disposable Prophy Angles (200 ct)

Request a demo at [nuprofreedom.com](http://nuprofreedom.com)



Save your patients and yourself from the risk of infection.  
Visit [dentsplysirona.com/preventive](http://dentsplysirona.com/preventive) today.



Dentsply Sirona Preventive offers a full range of solutions & support throughout the entire preventive procedure.

<sup>1</sup> Centers for Disease Control and Prevention. "Summary of Infection Prevention Practices in Dental Settings: Basic Expectations for Safe Care." Atlanta, GA: US Department of Health and Human Services, Centers for Disease Control and Prevention, National Center for Chronic Disease Prevention and Health Promotion, Division of Oral Health; March 2016.

<sup>2</sup> CDC MMWR Dec 19 2003. Guidelines for Infection Prevention Control in Dental Healthcare Settings. 2003.