

Long-term clinical evidence with Ankylos[®] Implant System

NEWLY PUBLISHED RESULTS WITH DENTSPLY SIRONA IMPLANTS PRODUCTS

12,737 Ankylos implants

4,206 Ankylos patients

followed up to

20 years

Ankylos implants show high implant survival rate and low rate of peri-implant bone loss after up to 20 years in function.

93.3%

Implant survival rate after 204 months

Unchanged implant design:

- Progressive threads
- Platform shift
- Morse taper connection



University study*



<1 mm**

* 36 oral surgeons treated patients at Frankfurt University, Germany

** Less than 1 mm vertical bone loss was observed for 85.2% of the implants after 17 years



Krebs M, Schmenger K, Neumann K, et al. Long-term evaluation of Ankylos dental implants, part i: 20-year life table analysis of a longitudinal study of more than 12,500 implants. Clin Implant Dent Relat Res 2015;17 Suppl 1:e275-86. [Abstract](#)



www.dentsplysirona.com

 **Dentsply
Sirona**
Implants