

# Atlantis® IO FLO compatibility chart

The Atlantis IO FLO is an impression component used during intraoral scanning to enable ordering of Atlantis products. Refer to instructions for use for further information, 4219-IFU Atlantis IO FLO\*.









Manufacturer	Implant	IO FLO code	REF
<b>BioHorizons</b>	Internal, Tapered Internal 3.5	IO-B-01	35212
	Internal, Tapered Internal 4.5	IO-B-02	35213
	Internal, Tapered Internal 5.7	IO-B-03	35214
<b>Biomet 3i</b>	Certain™ 3.25, 4/3	IO-F-01	35225
	Certain™ 4.0, 5.0, 6.0, 5/4	IO-F-02	35226
	Certain™ Prevail™ 3/4/3 (3.4), 4/3	IO-F-01	35225
	Certain™ Prevail™ 4/5/4, 5/4, 5/6/5, 6/5	IO-F-02	35226
	Certain™ XP 4/5, 5/6	IO-F-02	35226
<b>Camlog</b>	CONOLOG 3.3	IO-Q-01	35267
	CONOLOG 3.8	IO-Q-02	35268
	CONOLOG 4.3	IO-Q-02	35268
	CONOLOG 5.0	IO-Q-03	35269
	Screw-line Implant 3.3	IO-N-01	35239
	Screw-line/Root-line Implant 3.8	IO-N-02	35240
	Screw-line/Root-line Implant 4.3	IO-N-03	35241
	Screw-line/Root-line Implant 5.0 / 6.0	IO-N-04	35242
<b>Dentsply Sirona</b>	Ankylos® C/X 3.5, 4.5, 5.5, 7.0	IO-L-01*	35235
	Astra Tech Implant EV 3.0 <span style="color: green;">3.0 GREEN</span>	IO-P-01	35243
	Astra Tech Implant EV 3.6 <span style="color: purple;">3.6 PURPLE</span>	IO-P-02	35244
	Astra Tech Implant EV 4.2 <span style="color: yellow;">4.2 YELLOW</span>	IO-P-03	35245
	Astra Tech Implant EV 4.8 <span style="color: blue;">4.8 BLUE</span>	IO-P-04	35246
	Astra Tech Implant EV 5.4 <span style="color: brown;">5.4 BROWN</span>	IO-P-05	35247
	Astra Tech Implant EV Profile 4.2 <span style="color: yellow;">4.2 YELLOW</span>	IO-P-06	35248
	Astra Tech Implant EV Profile 4.8 <span style="color: blue;">4.8 BLUE</span>	IO-P-07	35249
	OsseoSpeed® / OsseoSpeed® TX 3.0 <span style="color: yellow;">3.0 YELLOW</span>	IO-C-01	35215
	OsseoSpeed® / OsseoSpeed® TX 3.5, 4.0 <span style="color: aqua;">3.5, 4.0 AQUA</span>	IO-C-02	35216
	OsseoSpeed® / OsseoSpeed® TX 4.5, 5.0 <span style="color: purple;">4.5, 5.0 LILAC</span>	IO-C-03	35217
	OsseoSpeed® TX Profile 4.5, 5.0	IO-C-04	35218
	Xive® S 3.0 <span style="color: red;">3.0</span>	IO-J-01	35230
	Xive® S 3.4 <span style="color: grey;">3.4</span>	IO-J-02	35231
	Xive® S 3.8 <span style="color: orange;">3.8</span>	IO-J-03	35232
	Xive® S 4.5 <span style="color: blue;">4.5</span>	IO-J-04	35264
	Xive® S 5.5 <span style="color: red;">5.5</span>	IO-J-05	35265

\*The Atlantis IO FLO for Ankylos C/X is for **single use only**. Refer to instructions for use for further information, 4190-IFU Atlantis IO FLO for Ankylos.

All trademarks and company names are the property of their respective owners.

Please turn over

# Atlantis® IO FLO Singles

Manufacturer	Implant	IO FLO code	REF
<b>MIS</b>	C1 NP 3.3	IO-S-02	34689
	C1 SP 3.75, 4.2	IO-S-03	34690
	C1 WP 5.0	IO-S-04	34694
	M4 NP 3.3	IO-B-04	35273
	M4 SP 3.75, 4.2	IO-B-01	35212
	M4 WP 5.0, 6.0	IO-B-02	35213
	SEVEN NP 3.3	IO-B-04	35273
	SEVEN SP 3.75, 4.2	IO-B-01	35212
	SEVEN WP 5.0, 6.0	IO-B-02	35213
	V3 NP 3.3	IO-S-01	34688
	V3 SP 3.9, 4.3, 5.0	IO-S-03	34690
<b>Nobel Biocare</b>	NobelActive™ 3.0	IO-K-03	35266
	NobelActive™ / NobelReplace™ Conical C NP 3.5 	IO-K-01	35233
	NobelActive™ / NobelReplace™ Conical C RP 4.3 / 5.0 	IO-K-02	35234
	NobelActive WP 5.5	IO-K-04	35272
	NobelReplace™ NP 3.5  (Tri-channel connection)	IO-D-01	35219
	NobelReplace™ RP 4.3  (Tri-channel connection)	IO-D-02	35220
	NobelReplace™ WP 5.0  (Tri-channel connection)	IO-D-03	35221
	NobelReplace™ 6.0  (Tri-channel connection)	IO-D-04	35222
<b>Osstem</b>	Osstem TS Mini 3.5	IO-R-01	35270
	Osstem TS Regular 4.0, 4.5, 5.0, 6.0, 7.0	IO-R-02	35271
<b>Straumann</b>	Bone Level 3.3 NC	IO-H-01	35228
	Bone Level 4.1 RC / 4.8 RC	IO-H-02	35229
<b>Zimmer Dental</b>	Screw-Vent™ / Screw-Vent™ friction fit 3.3, 3.7	IO-B-01	35212
	Screw-Vent™ / Screw-Vent™ friction fit 4.5	IO-B-02	35213
	Tapered Screw-Vent™ / Tapered Screw-Vent™ friction fit 3.5	IO-B-01	35212
	Tapered Screw-Vent™ / Tapered Screw-Vent™ friction fit 4.5	IO-B-02	35213
	Tapered Screw-Vent™ / Tapered Screw-Vent™ friction fit 5.7	IO-B-03	35214

All trademarks and company names are the property of their respective owners.