

Powder Opaque	Contact us regarding SKU information for Powder Opaque items	
Paste Opaque	3mL	
Paste Opaque - Bleach BL1/BL2	5360901220	
Paste Opaque - Bleach BL3/BL4	5360901230	
Paste Opaque - A1	5360903000	
Paste Opaque - A2	5360903010	
Paste Opaque - A3	5360903020	
Paste Opaque - A3.5	5360903030	
Paste Opaque - A4	5360903040	
Paste Opaque - B1	5360903050	
Paste Opaque - B2	5360903060	
Paste Opaque - B3	5360903070	
Paste Opaque - B4	5360903080	
Paste Opaque - C1	5360903090	
Paste Opaque - C2	5360903100	
Paste Opaque - C3	5360903110	
Paste Opaque - C4	5360903120	
Paste Opaque - D2	5360903130	
Paste Opaque - D3	5360903140	
Paste Opaque - D4	5360903150	
Opaque Modifier	3mL	
Paste Opaque - Neutral	5360903160	
Paste Opaque - Gum	5360903200	
Paste Opaque - Intensive Orange	5360903560	
Dentin	20g	75g
Dentin - Bleach	5360911203	5360911205
Dentin - Bleach BL1	5360911223	
Dentin - Bleach BL2	5360911233	
Dentin - Bleach BL3	5360911243	
Dentin - Bleach BL4	5360911253	
Dentin - A1	5360911003	5360911005
Dentin - A2	5360911013	5360911015
Dentin - A3	5360911023	5360911025
Dentin - A3.5	5360911033	5360911035
Dentin - A4	5360911043	5360911045
Dentin - B1	5360911053	5360911055
Dentin - B2	5360911063	5360911065
Dentin - B3	5360911073	5360911075
Dentin - B4	5360911083	5360911085
Dentin - C1	5360911093	5360911095
Dentin - C2	5360911103	5360911105
Dentin - C3	5360911113	5360911115
Dentin - C4	5360911123	5360911125
Dentin - D2	5360911133	5360911135
Dentin - D3	5360911143	5360911145
Dentin - D4	5360911153	5360911155
Action-i Dentine	20g	
Butter	5360911903	
Chocolate	5360911913	
Corn	5360911923	
Crème	5360911933	
Honey	5360911943	
Lavender	5360911953	
Marble	5360911963	
Effect Powders	20g	75g
Transpa Red	5360915103	
Transpa Clear	5360923003	5360923005
Gray Inside	5360915113	
Flu Inside 1	5360917003	
Flu Inside 2	5360917053	
White Surface, WS	5360923013	
Stand By, S By	5360923023	5360923025
Dentin Gingival	20g	
Gum 1	5360919003	
Gum 2	5360919013	
Gum 3	5360919023	
Gum 4	5360919033	
Power Chroma	20g	75g
Power Chroma - PC 1	5360955003	
Power Chroma - PC 2	5360955013	5360955015
Power Chroma - PC 3	5360955015	5360955025
Power Chroma - PC 4	5360955033	5360955035
Power Chroma - PC 5	5360955043	5360955045
Power Chroma - PC 6	5360955053	

Incisal	20g	75g
Incisal - Bleach BL1/BL2	5360921623	
Incisal - Bleach BL3/BL4	5360921633	
Incisal - S01	5360921503	5360921505
Incisal - S02	5360921513	5360921515
Incisal - S03	5360921523	5360921525
Incisal - S04	5360921533	5360921535
Incisal - S05	5360921543	5360921545
Incisal - S06	5360921553	5360921555
Opal Incisal	20g	75g
Opal Incisal - 1	5360927503	5360927505
Opal Incisal - 2	5360927513	5360927515
Opal Effect	20g	75g
Opal Effect - Sunrise	5360927003	5360927005
Opal Effect - Sunset	5360927023	5360927025
Opal Effect - Sky	5360927033	5360927035
Opal Effect - Ocean	5360927043	5360927045
Opal Effect - Fog	5360927053	5360927055
Add-On Corrector	20g	75g
Add-On Corrector - Final Kiss	5360937003	5360937005
Shoulder	20g	
Shoulder - Bleach	5360941223	
Shoulder - SM1	5360941003	
Shoulder - SM2	5360941013	
Shoulder - SM3	5360941023	
Shoulder - SM4	5360941033	
Shoulder - SM5	5360941043	
Final Shoulder	20g	
Final Shoulder - FSM1	5360941053	
Final Shoulder - FSM2	5360941063	
Final Shoulder - FSM3	5360941073	
Final Shoulder - FSM4	5360941083	
Final Shoulder - FSM5	5360941093	
Kit Configurations		
Full Assortment w/Paste Opaque (A-D shades)	EA	5360990131
Starter Set (A2, A3, A3.5, B1, B3, D3 shades)	EA	5360990141
Test Set (A2 shade)	EA	5360990151
Bleach Shades Set (BL1-4)	EA	5360990191
Artist Kit (Action-i Dentin)	EA	5360990161
Accessories		75g
Universal NP-Bonder, Paste	4mL	5368900501
NE-Bonder, Powder	20g	5360900503 5360900505
Ducera Liquid SD	50mL	5368272054
Ducera Liquid FORM	50mL	5368272034
Dentsply Sirona Stain and Glaze Liquid	15mL	601315
Kiss Shade Wheel	EA	22030
Dentsply Sirona Portioner Tool	EA	5369002400
Dentsply Sirona Universal Stain & Overglaze		
DS Universal Stain Kit (all items below)	EA	600700
DS Universal Body Stain - S0	5g	605520
DS Universal Body Stain - S1	5g	605521
DS Universal Body Stain - S2	5g	605522
DS Universal Body Stain - S3	5g	605523
DS Universal Body Stain - S4	5g	605524
DS Universal Incisal Stain - i1	5g	605531
DS Universal Incisal Stain - i2	5g	605532
DS Universal Stain - White	5g	605500
DS Universal Stain - Crème	5g	605501
DS Universal Stain - Sunset	5g	605502
DS Universal Stain - Copper	5g	605503
DS Universal Stain - Khaki	5g	605504
DS Universal Stain - Olive	5g	605505
DS Universal Stain - Mahogany	5g	605506
DS Universal Stain - Violet	5g	605507
DS Universal Stain - Raspberry	5g	605508
DS Universal Stain - Purple	5g	605509
DS Universal Stain - Grey	5g	605510
DS Universal Stain - Chestnut	5g	605511
DS Universal Stain - Blue	5g	605512
DS Universal Stain - Pink	5g	605513
DS Universal Overglaze	5g	605540
DS Universal Overglaze - High Flu	5g	605542

DequDent GmbH
 Rodenbacher Chaussee 4
 63457 Hanau-Wolfgang
 Germany
 +49 6181 59-50
 dentsplysirona.com/ceramics

Dentsply Sirona Prosthetics
 570 West College Avenue
 York, PA 17401
 800-243-1942

©2018 Dentsply Sirona. All rights reserved.

Form No. 22319 Rev. 2018-06



Duceram® Kiss
 Metal Veneering Porcelain

Craft meets creativity



Duceram® Kiss (Kiss System)

Craft meets creativity

Whether Duceram Kiss, Duceragold® or Cercon® ceram Kiss - the complete Kiss system offers simplicity and efficiency across framework materials.

- + Multi-framework concept
- + Perfect aesthetics
- + Simple and secure

Duceram Kiss is a high-fusing ceramic material for veneering crowns and bridges with frameworks made of dental alloys with a CTE range of 13.8 to 15.4 $\mu\text{m}/\text{m} \cdot \text{K}$ (25 - 600 °C).


Duceram Kiss is suitable for veneering precious and non-precious metal frameworks with a coefficient of thermal expansion (CTE) from 13.8 to 15.4 $\mu\text{m}/\text{m} \cdot \text{K}$ (25 - 600 °C).



KISS SHADE WHEEL

Shade accurate material selection made easy

The Kiss Shade Wheel lets you quickly and easily identify the respective Kiss ceramic materials corresponding to all 16 VITA A1-D4 shades.*



*VITA and VITA classical A1-D4 are registered trademarks of VITA Zahnfabrik H. Rauter GmbH & Co. KG


Laboratory work: Joachim Wille

Multi-framework efficiency

All Kiss systems share the same number of components for a consistent veneering concept - just 73 masses in total. Four Bleach (BL1-BL4) and 16 VITA classical dentin shades and associated aesthetic masses are included in the Kiss system.

ARTIST KIT

Looking for more advanced veneering possibilities?



Duceram Kiss Artist Kit introduces additional masses for even greater individualization supporting more advanced veneering concepts with rich colors including new Gum shades.

UNIVERSAL FINISHING SYSTEM

For metal-based and all-ceramic restorations



Dentsply Sirona Universal Stains' broad compatibility eliminates multiple finishing systems in your lab - learn more at dentsplysirona.com/ceramics



Duceram Kiss veneering porcelain

Duceram Kiss is the premium, high-fusing veneering system for a broad range of classic metal frameworks, including precious and non-precious alloys. The most distinctive characteristics of Duceram Kiss are its exceptional shade fidelity, robust and simple handling and firing stable opal incisals, ensuring perfect aesthetics.



Multi-framework concept

+ Common veneering concept across multiple frameworks creates work efficiency using just 73 masses saving time and money.

Perfect aesthetics

+ Optimized pigments enhance shade fidelity and a patented process maintains exceptional opalescent effect in a high-fusing ceramic.

Simple and secure

+ Secure process utilizes just one simple layering system for basic and highly individualized techniques across multiple framework options.