



Precise Replication... Down To The Last Detail.

Trixa is the putty that enables you to replicate every detail of your removable appliance from festooning and stippling to exact tooth location.

Trixa is ideal for:

- Repairs, relines, rebases, duplicating dentures and cold curing partials
- Various matrices or positioning guides
- Insulating putty use when flasking
- Removable gingival masks (for implants and crown and bridge)
- Indexing teeth

Trixa Times	
Task	Time
Mixing	45 seconds
Working	1-2 minutes
Setting	4-7 minutes