

Palodent® 360

Circumferential Matrix System

Simplicity comes full circle

Palodent® 360 is a circumferential matrix with integrated tightening mechanism and ultra-thin, high strength matrix band - no need for bulky retainers or applicators. The contoured matrix provides anatomical shape and optimal contacts with a simple twist of your fingers.

Faster and easier restorations than you ever thought possible thanks to the integrated tightener/retainer which allows for:

- Efficiency (no set up required)
- Better access
- Clear view of the operating field
- Greater patient comfort
- Ease of placement




Placement versatility


Palodent® 360 system is designed for challenging Class II cases where a sectional matrix cannot be used and for single-unit crown build ups. The tightening wheel can be positioned where most convenient; lingual, buccal or a more proximal aspect.

Examples of clinical indications:


Missing cusp




Missing adjacent tooth



Misaligned tooth




Tooth isolation for core build-up



Palodent® 360

Natural contours and optimal contacts with a simple twist of the fingers

The innovative design means a twist of your fingers is all it takes to achieve tight contacts with anatomically correct contours.



Traditional matrix band

- ✗ Bulky system
- ✗ Failure to restore proximal anatomy, increasing the risk of fracture, recurrent caries and periodontal disease
- ✗ Weak or open contact point - high risk of food impaction
- ✗ Extensive finishing required

Palodent® matrix system restoration

- ✓ Easy to use
- ✓ Natural contours
- ✓ Tight, anatomically correct contacts
- ✓ Minimal interproximal finishing

Palodent® Matrix Systems

659600v Complete System

- 100 Palodent® V3 Matrices – 25 each size: 3.5mm, 4.5mm, 5.5mm, 6.5mm
- 75 Palodent® V3 Wedges – 25 each size: Small, Medium, Large
- 30 Palodent® V3 WedgeGuards – 10 each size: Small, Medium, Large
- 1 Palodent® V3 Universal Ring
- 1 Palodent® V3 Narrow Ring
- 1 Palodent® V3 Forceps
- 1 Palodent® V3 Pin Tweezers
- 24 Palodent® 360 Matrices – 5.5mm



Palodent® 360

Circumferential Matrix System

Matrix Band Refill (48)

- 659601 4.5mm – Dark Blue
- 659602 5.5mm – Medium Blue
- 659603 6.5mm – Light Blue

Wedge Refill (100)

- 659780v Small – Dark Blue
- 659790v Medium – Medium Blue
- 659800v Large – Light Blue

WedgeGuard Refill (50)

- 659830v Small – Dark Blue
- 659840v Medium – Medium Blue
- 659850v Large – Light Blue

www.dentsplysirona.com