

# X-SMART<sup>®</sup>

## ENDO MOTOR

**Compact.**

**Portable.**

**Easy to use.**



**DENSPLY**  
TULSA DENTAL  
SPECIALTIES

# X-SMART®

## ENDO MOTOR

**Compact. Portable. Easy to use.**

### Technical Features:

- Operates on battery or electrical power
- Battery recharging time: approximately 5 hours
- Time of use on battery: approximately 2 hours
- Supplied with a 16:1 compact contra angle
- Operates with or without a pedal (sold separately)
- Available alone or with file and obturation starter systems



Hand piece with integrated on/off switch  
(only 3.22 ounces or 92 grams)

Large LCD screen

Speed range\*:  
120 to 800 rpm

Auto reverse mode

Torque control range\*:  
0.6 to 5.2 N-cm  
(61 - 102 g-cm)

9 adjustable  
program selections

\*with a supplied 16:1  
compact contra angle



**Contact your Tulsa Dental Specialties representative for more information.**

1-800-662-1202  
For the latest information consult  
[www.TulsaDentalSpecialties.com](http://www.TulsaDentalSpecialties.com)

**DENTSPLY**  
TULSA DENTAL  
SPECIALTIES

DESCRIPTION	SKU
X-Smart® Endo Motor (Complete Kit)	XSMART
X-Smart® Contra Angle	XSMARTCA
X-Smart® Foot Pedal	XSMARTFP
X-Smart® Motor w/ Complete ProSystem GT® Rotary & Obturation System	PROSYSGTXSMT
X-Smart® Motor w/ ProFile® GT® Rotary System	XSMTSYSGT
X-Smart® Motor w/ Complete ProTaper® Universal Rotary & Obturation System	PTUENDOSYSXSM
X-Smart® Motor w/ ProTaper® Universal Rotary System	XSMARTPTU
X-Smart® Motor w/ 4 Packs of ProTaper® Universal Files	PTUSYSXSMT