

Compatibility list Flex/Green/White Lid adapters for DAC Universal GUI/Text/D/S

Manufacturer	Current portfolio	Instrument	Ref. manufacturer	Type	Indictaion	Adapter Ref. Sirona	Adapter Sirona US Scaler tip	Adapter Sirona US Scaler handpiece (Sirosonic / L)	Adapter Sirona US Scaler handpiece (Sirosonic TL)	Adapter EMS US Scaler tip (Piezon®, Piezon® LED)	Adapter EMS US Scaler handpiece (Piezon®, Piezon® LED)	Adapter KaVo SONICflex Handpiece (2003/L, 2008/L)	Adapter KaVo SONICflex 2003 tip adapter (only for SONICflex 2003/L tips)	Adapter Satelec US Scaler tip	Adapter Sirona Sprayvit nozzle G	Adapter Sirona Sprayvit nozzle 4000	Adapter EMS AIR-FLOW-HANDY 3.0	Adapter Air-Flow Handpiece	Adapter Satelec Slim	Adapter Satelec Newtron	Adapter Satelec Newtron LED	Adapter X-Smart Pro/Pro+ sleeve
EMS	yes	Instrument A	DS-001 A	Tip	Scaling	6610708				x												
EMS	yes	Instrument P	DS-011 A	Tip	Scaling	6610708				x												
EMS	yes	Instrument PS	DS-016 A	Tip	Scaling	6610708				x												
EMS	yes	Instrument B	DS-003 A	Tip	Scaling	6610708				x												
EMS	yes	Instrument C	DS-002 A	Tip	Scaling	6610708				x												
EMS	yes	Instrument PL1	DS-030 A	Tip	Perio	6610708				x												
EMS	yes	Instrument PL2	DS-031 A	Tip	Perio	6610708				x												
EMS	yes	Instrument PL3	DS-032 A	Tip	Perio	6610708				x												
EMS	yes	Instrument PL4	DS-033 A	Tip	Perio	6610708				x												
EMS	yes	Instrument PL5	DS-034 A	Tip	Perio	6610708				x												
EMS	yes	Instrument HPL3	DS-075 A	Tip	Perio	6610708				x												
EMS	yes	Instrument DPL3	DS-065 A	Tip	Perio	6610708				x												
EMS	yes	Piezon® LED Handpiece	EN-060/A	Handpiece	n.a.	6613538					x											
EMS	yes	Piezon® Handpiece	EN-061/A	Handpiece	n.a.	6613538					x											
EMS	yes	Piezon® Handpiece	EN-041	Handpiece	n. a.	6750090					x											
EMS	yes	Piezon® Handpiece	EN-046	Handpiece	n. a.	6750090					x											
EMS	yes	AIRFLOW MAX Handpiece	EL-308/XC/YC	Handpiece	n.a.	6820356												x				
EMS	yes	AIRFLOW Handpiece	EL-308/LC/C	Handpiece	n.a.	6820356												x				
EMS	yes	PERIOFLOW Handpiece	EL-354	Handpiece	n.a.	6820356												x				
EMS	yes	AIRFLOW PLUS Handy 3.0 Perio	EL-545/A	Handpiece nozzle	n.a.	6623453											x					
EMS	yes	AIRFLOW PERIO Handy 3.0	EL-542	Handpiece nozzle	n.a.	6623453											x					
EMS	yes	AIRFLOW Handy 3.0	EL-540	Handpiece nozzle	n.a.	6623453											x					
KaVo	yes	KaVo PiezoLED	1.007.3995	Handpiece	n.a.	6613538					x											
KaVo	yes	KaVo SONICflex 2003**	1.000.4246	Handpiece	n.a.	6732056						x										
KaVo	yes	KaVo SONICflex LUX 2003 L**	1.000.4243	Handpiece	n.a.	6732056						x										
KaVo	yes	KaVo SONICflex quick 2008**	1.005.9311	Handpiece	n.a.	6732056						x										
KaVo	yes	KaVo SONICflex quick 2008 L**	1.005.9310	Handpiece	n.a.	6732056						x										
KaVo	yes	KaVo SONICflex 2003 Spitzen**	several	Tip	check Dental-Information 09/2020	6735646							x									
Sirona	yes	Ultrasonic tip 1L	8924250	Tip	Scaling	6536127	x															
Sirona	yes	Ultrasonic tip 2L	8924268	Tip	Scaling	6536127	x															
Sirona	yes	Ultrasonic tip 3L	1894455	Tip	Scaling	6536127	x															
Sirona	yes	Ultrasonic tip 4L	1894463	Tip	Scaling	6536127	x															
Sirona	yes	Ultrasonic tip SiroRetro KL	5457887	Tip	Endodontics	6536127	x															
Sirona	yes	Ultrasonic tip SiroRetro KR	5457895	Tip	Endodontics	6536127	x															
Sirona	yes	Ultrasonic tip SiroRetro GL	5457903	Tip	Endodontics	6536127	x															
Sirona	yes	Ultrasonic tip SiroRetro GR	5457911	Tip	Endodontics	6536127	x															

