

Cavitron® Scaling

Enhance your scaling experience with the latest innovations in Cavitron® 30k Ultrasonic Scaling Systems.

Cavitron® 30k units offer features that help provide you and your patients the best possible dental care through innovations in ergonomics, infection prevention and insert selection.

30,000
strokes per second

VS.


25,000
strokes per second



Infection Prevention

Handpiece

Cavitron Steri-Mate® handpieces are detachable and autoclavable providing assurance of control against cross-contamination.



Handpiece

The 25K Bobcat handpiece is not detachable or autoclavable.



Ultrasonic scaling handpieces are considered semi-critical and should be heat sterilized between patients to reduce cross-contamination.

Ergonomics

Wireless Foot Pedal*

Wireless 360° foot pedal* is ergonomically designed to provide ease of use and mobility.



Wired Foot Pedal

The 25K Bobcat has a square, wired foot pedal.

Rotating Handpiece**

Cavitron Steri-Mate® 360** fully rotating handpiece allows for free flowing movement of an ultrasonic insert within the oral cavity.



Nearly one-third of dental health professionals cite musculoskeletal disorders as the reason for early retirement.¹ Ergonomic features ensure an intuitive performance and allow clinicians to focus on patient care.

Handpiece

The 25K Bobcat handpiece is longer and heavier.



Innovation

Wide Range of Inserts with Grip Options

- General Scaling
- Implant Maintenance
- Periodontal Debridement
- Surgical Applications



Just like with hand scaling, inserts are built for removing different levels of deposit and specific types of tooth anatomy. Selecting and using the proper inserts can help improve patient comfort and efficiency of your practice.

Limited Inserts

Limited tip design and finger grip options.



1. OSHA Watch Updates July/August 2007 Volume 9, No. 4 page 2 – References The Quality America Dental OSHA Safety Program Manual Tab 4
 * Available on Cavitron Plus and 300 models
 ** Compatible ONLY with Cavitron 300 Series and Integrated G139 Systems

Ordering Information

Cavitron® 300 Series

- 8161429 Cavitron® 300 Series Ultrasonic Scaling System Package
 - (1) Cavitron® 300 Series Ultrasonic Scaling System Unit
 - (3) Cavitron® Steri-Mate® 360 Handpieces
 - (1) Cavitron® Slimline 10S Fitgrip® 30K Ultrasonic Insert
 - (3) Cavitron® Inserts of your choice*
 - 3 Years of Cavitron Care: Maintenance Plan
 - 4-Year Warranty
- 8270310 Cavitron® 300 Series Ultrasonic Scaling System Unit
 - (1) Cavitron® 300 Series Ultrasonic Scaling System Unit
 - (1) Cavitron® Steri-Mate® 360 Handpiece
 - 3-Year Warranty

Cavitron® Jet Plus Ultrasonic Scaler and Air Polishing System

- 8187501 Cavitron® Jet Plus Ultrasonic Scaler and Air Polishing System
- 8161428 Cavitron® Jet® Plus Package

Cavitron® Plus Ultrasonic Scaling System

- 8184001 Cavitron® Plus Ultrasonic Scaling System
- 8161425 Cavitron® Plus Package

Cavitron® Prophy-Jet Air Polishing System

- 8195001 Cavitron® Prophy Jet® Air Polishing System
- 8161427 Cavitron® Prophy Jet® Package

Cavitron® Select® SPS™ Ultrasonic Scaling System

- 81302 Cavitron® Select® SPS™ Ultrasonic Scaling System
- 81308 Cavitron® Select® SPS™ Ultrasonic Scaling System with Reservoir
- 8161424 Cavitron® Select® SPS™ Package

Cavitron® Slimline® Ultrasonic Inserts

- 80395 Cavitron® Slimline® 10S 30K Ultrasonic Insert
- 82005 Cavitron® Slimline® 10S Fitgrip® 30K Ultrasonic Insert
- 80397 Cavitron® Slimline® 10R 30K Ultrasonic Insert
- 82007 Cavitron® Slimline® 10R Fitgrip® 30K Ultrasonic Insert
- 80396 Cavitron® Slimline® 10L 30K Ultrasonic Insert
- 82006 Cavitron® Slimline® 10L Fitgrip® 30K Ultrasonic Insert
- 81570 Cavitron® Slimline® 1000 30K Ultrasonic Insert
- 82008 Cavitron® Slimline® 1000 Fitgrip® 30K Ultrasonic Insert

Cavitron® Powerline® Ultrasonic Inserts

- 80797 Cavitron® Powerline® 3 30K Ultrasonic Insert
- 82001 Cavitron® Powerline® 3 Fitgrip® 30K Ultrasonic Insert
- 80294 Cavitron® Powerline® 10 30K Ultrasonic Insert
- 82002 Cavitron® Powerline® 10 Fitgrip® 30K Ultrasonic Insert
- 80798 Cavitron® Powerline® 100 30K Ultrasonic Insert
- 82003 Cavitron® Powerline® 100 Fitgrip® 30K Ultrasonic Insert
- 80799 Cavitron® Powerline® 1000 30K Ultrasonic Insert
- 82004 Cavitron® Powerline® 1000 Fitgrip® 30K Ultrasonic Insert

Cavitron® Thinsert® Ultrasonic Inserts

- 81551 Cavitron® Thinsert® 30K Ultrasonic Insert
- 82009 Cavitron® Thinsert® Fitgrip® 30K Ultrasonic Insert

Cavitron® Steri-Mate® 360 Rotating Handpiece

- 8183201 Cavitron® Steri-Mate® 360 Handpiece - Single Pack
- 8183203 Cavitron® Steri-Mate® 360 Handpiece - 3-Pack

For more information or to request a demo, call **1.800.989.8826** or visit dentsplysirona.com/preventive.

Cavitron® Inserts

For hygienists who strive for efficiency, comfort, and performance, the Cavitron® family of inserts work together to deliver the highest standard in ultrasonic scaling.

A comprehensive line of ultrasonic scaling inserts designed to best serve you, your patients, and your practice.



Thinsert®

Thinnest tip for maximum access in difficult areas

Slimline®


Slim design for improved subgingival access and accessibility

Powerline®

Standard size for efficient removal of heavier deposits




Set up



Diagnosis



Anesthetics



Scaling



Polishing



Treatments

Dentsply Sirona Preventive offers a full range of solutions & support throughout the entire preventive procedure.