

Cavitron®

Tabletop System Comparison

Scaling

Air Polishing



Cavitron® Select™ SPS™



Cavitron® Plus



Cavitron® 300 Series



Cavitron® Prophy Jet®



Cavitron® Jet Plus

	Portable & Compact	Improved comfort & efficiency	The ultimate scaling experience¹	Air Polishing	Combo Scaling & Air Polishing
Scaling and Air Polishing Capabilities	●	●	●	●	✓
Prophy Mode Auto Cycles	●	●	●	✓	✓
Wireless Tap-On Technology Foot Pedal	●	✓	✓	✓	✓
Rinse Function	●	✓	✓	✓	✓
Purge Mode	●	✓	✓	✓	✓
Pivoting Handpiece Cable	✓	✓	✓	✓	✓
Water Control on Handpiece	✓	✓	✓	✓	✓
Boost Mode	✓	✓	✓	●	✓
Intelligent Feedback Control (SPS Technology)	✓	✓	✓	●	✓
Steri-Mate/Jet-Mate Autoclavable Detachable Handpiece	✓	✓	✓	✓	✓
Steri-Mate 360 Autoclavable Detachable Swivel Handpiece	●	●	✓	●	●
Touch Screen Interface	●	●	✓	●	●
Memory Presets	●	●	✓	●	●
Featherweight Detachable Handpiece Cord	●	●	✓	●	●
Water Filter Indication Light	●	●	✓	●	●
Optional Reservoir Pump/Bottle System (NOTE: SKU# 90146 is sold separately and is only compatible with Cavitron Select SPS)	SKU# 90146	●	●	●	●
Warranty	1 Year	2 Years	3 Years	2 Years	2 Years
Dentsply Sirona Part Number	81302	8184001	8270310	8195001	8187501
MSRP					
Dentsply Sirona Rebate					
Net MSRP					

Package Options

	Select™ SPS™ Package	Plus Package	300 Series Package	Prophy Jet® Package	Jet Plus Package
Warranty	2 Years	3 Years	4 Years	3 Years	3 Years
Cavitron Care Maintenance Plan	Not Available	2 Years	3 Years	2 Years	3 Years
Additional Handpieces	●	●	(2) Steri-Mate 360 Handpieces	●	(2) Jet-Mate Handpieces
Cavitron Inserts (NOTE: (1) SlimLine 10S FitGrip included in package, Redeem 3 of your choice. Not applicable to the Prophy Jet Package)	(4) Cavitron 30K Ultrasonic Inserts	(4) Cavitron 30K Ultrasonic Inserts	(4) Cavitron 30K Ultrasonic Inserts	(2) Cavitron Air Polishing Inserts	(4) Cavitron 30K Ultrasonic Inserts
Dentsply Sirona Part Number	8161424	8161425	8161429	8161427	8161428
Package MSRP					
Dentsply Sirona Rebate					
Net Package Price					
Product Retail Value (Inserts and/or Handpieces)					
Cavitron Care Value (3 Year Warranty and/or Maintenance Plan)	\$150	\$550	\$750	\$550	\$750
Total Package Value (Rebate + Product Value + Cavitron Care Value)					

Cavitron Care Maintenance Plan

Available exclusively with Cavitron Packages!

Cavitron Care will send you a maintenance reminder annually to schedule service on your Cavitron for the following items:

Unit Model	Scheduled Service Included
Cavitron 300 Series (G310)	3 years
Cavitron Jet Plus (G137)	3 years
Cavitron Plus (G136)	2 years
Cavitron Prophy Jet (G138)	2 years
Cavitron Select SPS (G124)	Not Available

General Operation and Function:	Electrical System:	Water Systems:	Air System: (Jet Plus & Prophy Jet Systems)
<ul style="list-style-type: none"> Evaluation and verification of all clinical functions Evaluation and repair of malfunctions Installation of recommended factory upgrades 	<ul style="list-style-type: none"> Calibration of ultrasonic circuits Test, verify and repair all electrical functions as needed Inspection and repair of electrical cable connections 	<ul style="list-style-type: none"> Inspection, repair/replacement of water control valves Testing and calibration of water flow with replacement of clogged or damaged components Evaluation and/or replacement of water supply line filter Repair or replacement of any components found to be leaking 	<ul style="list-style-type: none"> Evaluation and/or replacement of control valves Replacement of internal tubing and filters Testing and calibration of air flow rate with replacement of needed components Evaluation and/or replacement of main air supply line filter Repair or replacement of any components found to be leaking air

Cavitron® Ultrasonic Inserts

Each Cavitron Package includes (1) Slimline 10S FitGrip insert. Choose (3) additional inserts from the options below!

Slimline®



10 - Left Curve
Low to Medium Power
30K 80396
30K Fitgrip® 82006



10 - Right Curve
Low to Medium Power
30K 80397
30K Fitgrip® 82007



10 - Straight
Low to Medium Power
30K 80395
30K Fitgrip® 82005



1000
Low to High Power
30K 81570
30K Fitgrip® 82008

Powerline™



3
Low to High Power
30K 80797
30K Fitgrip® 82001



10
Low to High Power
30K 80294
30K Fitgrip® 82002



100
Low to High Power
30K 80798
30K Fitgrip® 82003



1000
Low to High Power
30K 80799
30K Fitgrip® 82004

When selecting inserts, consider the following:

Type of deposit

Inserts are uniquely made to handle different levels of deposit. Consider the force needed to remove the deposits.



Thicker inserts with larger tips are designed for heavier levels of calculus.



Thinner inserts with slim, adaptable tips are designed for the removal of light to moderate deposits or biofilm.

Treatment site

Select inserts that are specially designed for the type of tooth anatomy you will be working on.



Straight inserts are designed for flat tooth anatomy.



Curved inserts are designed to give you the best results in furcations and concavities.

Thinsert®



Thinsert
Low to High Power
30K 81551
30K Fitgrip® 82009

Simplified registration and redemption at: cavitronsystems.com