

Your individual options

Here you will find the complete list of all individual options and you can see which equipment is already part of the basic units.

Basic unit equipment options

VB CS TS

Options for the patient chair	VB	CS	TS
Double-articulated headrest	■	■	■
Motor-driven headrest	□	□	□
Premium upholstery	■	■	■
Premium thermo upholstery	□	□	□
Lounge upholstery	□	□	□
Lumbar support function	□	□	□
Swiveling armrest, right or left	□	□	□
4-way foot control	■	■	■
Foot control, cable version	■	■	■
Wireless foot control (Smart Control or Classic version)	□	□	□
External device connection	□	□	□
Booster cushion for small patients (C seat cushion)	□	□	□
Installation plate to level out uneven floors	□	□	□

Options for the dentist element (5 handpiece positions, removable handle & instrument holder)	VB	CS	TS
Sprayvit M syringe with light	□	□	□
1 st drive (for motor or turbine)	■	■	■
2 nd drive (for motor or turbine)	■	■	■
3 rd drive (for motor or turbine)	□	□	□
BL electric motor	□	□	□
BL ISO C electric motor	□	□	□
BL ISO E electric motor	□	□	□
BL Implant electric motor	□	□	□
Implant 20:1 contra-angle handpiece	□	□	□
Turbine hose	□	□	□
SiroSonic TL ultrasonic handpiece	□	□	□
Hands-free operation with cursor control	■	■	■
Implantology function	□	□	□
Saline solution pump	□	□	□
Endodontic treatment function	□	□	□
Contra-angle handpiece T1 Line Endo 6 L with light	□	□	□
Reciprocal file systems	□	□	□
ApexLocator function with display on the Sivation monitor	□	□	□
Swiveling tray holder for 1 or 2 standard size trays		□	□
Panoramic X-ray image viewer			□

Options for the assistant element (with 4 positions, removable holder, and silicone mat)	VB	CS	TS
Satelec mini LED curing light	□	□	□
Surgical suction device	□	□	□
2 nd HVE kit	□	□	□
Sprayvit M syringe with light	□	□	□
Suction device	■	■	■
Saliva ejector	■	■	■

Basic unit equipment options

VB CS TS

Options for the water unit	VB	CS	TS
Disinfection device with permanent disinfection and sanitization function	■	■	■
Chemical suction hose cleaning	■	■	■
Cuspidor valve for wet suction system	□	□	□
Multi-station valve for wet suction system	■	■	■
Automatic separator	□	□	□
Amalgam separator	□	□	□
Swiveling cuspidor	■	■	■
Motorized cuspidor swivel	□	□	□
Without cuspidor	□	□	□
Tray with clip-on tray and handle	□		

Options for Sivation digital	VB	CS	TS
SiroCam AF	□	□	□
SiroCam AF+	□	□	□
22" monitor on the tray incl. tray	□		
22" monitor on the light support column	□	□	□
Additional holder for camera	□		□
Camera supply in the dentist element	□	□	□
Camera supply in the assistant element		□	□
Sivation ready installation kit with USB interface	□	□	□
Cable set for PC connection	□	□	□

Working stools (optional long or short gas spring)	VB	CS	TS
Hugo Freehand (with foot release)	□	□	□
Hugo Manual (with manual release)	□	□	□
Hugo Manual Plus (with manual release and foot ring)	□	□	□
Carl Manual (with manual release)	□	□	□
Carl Manual Plus (optional long or short gas spring)	□	□	□
Theo (w/saddle-style seat)	□	□	□

Operating light/Heliident plus options	VB	CS	TS
LEDview Plus operating light (triple-jointed support arm, with sensor)	□	□	□
LEDview Plus ceiling model (with sensor)	□	□	□
LEDview Plus / Heliident Plus ceiling combination	□	□	□
Adapter for Heliident Plus unit model	□	□	□

VB Sinius (sliding track version)
 CS Sinius CS (whip arm version)
 TS Sinius TS (OTP version)
 ■ Included in the basic unit
 □ Optional